



FOR FURTHER INFORMATION CONTACT
USL Customer Services
FP 0800 658 814 FF 0800 830 660
e. customerservices@uslmedical.co.nz

# Biatain® Silicone Superior absorption - faster healing



The Biatain Silicone range (standard & Lite), delivers a secure fit without compromising the superior absorption of Biatain foam. Perforated, soft silicone adhesive wound contact layer ensures a secure fit to keep the dressing in place whilst delivering minimal pain / trauma upon removal. This silicone technology allows Biatain foam to absorb exudate and conform to the wound bed.

### Biatain® Silicone provides the following benefits:

- Superior absorption
- Secure fit
- · Non traumatic removal
- Ease of application
- · Latex free



Biatain® Silicone Shapes and Sizes

### Biatain® Silicone Adhesive Foam

Biatain Silicone is indicated for a wide range of exuding wounds when gentle adhesion is required.



| Code  | USL<br>Code | Size cm*                | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|-------|-------------|-------------------------|-------------------------|------------------|---------------------|
| 33434 | 13286       | 7.5 x 7.5               | Sterile                 | 10               | 320                 |
| 33435 | 13287       | 10 x 10                 | Sterile                 | 10               | 240                 |
| 33436 | 13288       | 12.5 x 12.5             | Sterile                 | 10               | 240                 |
| 33437 | 13289       | 15 x 15                 | Sterile                 | 5                | 120                 |
| 33438 | 13290       | 17.5 x 17.5             | Sterile                 | 5                | 80                  |
| 33400 | 14790       | 10 x 20                 | Sterile                 | 5                | 200                 |
| 33401 | 14791       | 10 x 30                 | Sterile                 | 5                | 100                 |
| 33404 | 14792       | 15 x 19 Sacral          | Sterile                 | 5                | 80                  |
| 33405 | 14793       | 25 x 25 Sacral          | Sterile                 | 5                | 90                  |
| 33406 | 14794       | 18 x 18 Heel            | Sterile                 | 5                | 40                  |
| 33408 | 14795       | 14 x 19.5<br>Multishape | Sterile                 | 5                | 80                  |

<sup>\*</sup>Indicates dressing size, for pad size refer to retail box

### Biatain® Silicone Lite Adhesive Foam

Biatain Silicone Lite is indicated for lightly exuding wounds when gentle adhesion is required. Its thinner composition offers more comfort and mobility to the patient.



| Code  | USL<br>Code | Size cm*    | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|-------|-------------|-------------|-------------------------|------------------|---------------------|
| 33444 | 33444       | 7.5 x 7.5   | Sterile                 | 10               | 320                 |
| 33445 | 33445       | 10 x 10     | Sterile                 | 10               | 240                 |
| 33446 | 33446       | 12.5 x 12.5 | Sterile                 | 10               | 240                 |
| 33452 | 14796       | 5 x 5       | Sterile                 | 5                | 120                 |
| 33453 | 14798       | 5 x 12.5    | Sterile                 | 5                | 140                 |

<sup>\*</sup>Indicates dressing size, for pad size refer to retail box

### Biatain<sup>®</sup>

### Superior absorption - faster healing

The Biatain range offers you a selection of foam dressings for a wide range of exuding wounds. All dressings are proven to deliver superior exudate management with minimal risk of leakage and maceration. It is well known that patient comfort is crucial for successful wound healing. Having listened to our customers, Biatain is now offering even more patient-friendly solutions.



#### 3D polymer structure

- Superior absorption and retention even under compression
- · Close wound conformability
- · Optimal moist wound environment

#### User friendly design

- Discreet appearance
- · Latex free

### Bevelled edges

- Significantly lower risk of pressure marks
- Fully sealed edges

The Biatain range provides a unique combination of very high patient comfort and superior exudate management across all dressing solutions.

### Biatain® Non-Adhesive Foam

Biatain non-adhesive is ideal for exuding wounds where the surrounding skin is fragile. Designed for superior comfort, this dressing combines all the benefits of the Biatain foam range to deliver an exceptional patient experience.



| Code | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|---------|-------------------------|------------------|---------------------|
| 6105 | 13276       | 5 x 7   | Sterile                 | 10               | 240                 |
| 3410 | 13277       | 10 x 10 | Sterile                 | 10               | 240                 |
| 3412 | 13278       | 10 x 20 | Sterile                 | 5                | 160                 |
| 3413 | 13279       | 15 x 15 | Sterile                 | 5                | 80                  |
| 3416 | 13280       | 20 x 20 | Sterile                 | 5                | 80                  |

### Biatain® Adhesive Foam

Biatain adhesive foam dressing has a reduced thickness, softer foam and discreet design with the same absorbency. The dressing uses a Hydrocolloid product as the adhesive component, ensuring skin friendliness.



| Code | USL<br>Code | Size cm*    | Sterile/<br>Non Sterile |    | Units Per<br>Carton |
|------|-------------|-------------|-------------------------|----|---------------------|
| 3462 | 13281       | 7.5 x 7.5   | Sterile                 | 10 | 320                 |
| 3420 | 13282       | 12.5 x 12.5 | Sterile                 | 10 | 240                 |
| 3423 | 13283       | 18 x 18     | Sterile                 | 5  | 80                  |

<sup>\*</sup>Indicates dressing size, for pad size refer to retail box

### Biatain® Adhesive Foam Shapes

The unique shapes allow the dressings to conform to the heel/sacrum. The dressing uses a Hydrocolloid product as the adhesive component, ensuring skin friendliness.

| and the same of the |
|---------------------|
|                     |
| and the second      |
|                     |


|   | Code | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|---|------|-------------|---------|-------------------------|------------------|---------------------|
|   | 3488 | 13285       | 19 x 20 | Sterile Heel            | 5                | 80                  |
| 3 | 3485 | 13284       | 23 x 23 | Sterile Sacral          | 5                | 80                  |

### Biatain® Soft-Hold Foam

Biatain Soft-Hold is an additional variant in our non-adherent range. It enables you to situate and keep the dressing in place - with a soft touch, while securing a secondary bandage and applying compression therapy - with minimal trauma to surrounding skin. Gentle adhesive loops of elastomer on the face of the foam ensure dressing security and turn into a soft gel on contact with exudate, creating an optimal moist wound healing environment and traumatic removal.




| Code | USL<br>Code | Size cm | Sterile/<br>Non Sterile |   | Units Per<br>Carton |
|------|-------------|---------|-------------------------|---|---------------------|
| 3470 | 13292       | 10 x 10 | Sterile                 | 5 | 240                 |
| 3472 | 13293       | 10 x 20 | Sterile                 | 5 | 160                 |
| 3475 | 13291       | 15 x 15 | Sterile                 | 5 | 80                  |



## Biatain® Alginate

Biatain Alginate is made of a unique combination of alginate and carboxymethylcellulose (CMC) for superior absorption and enhanced dressing integrity i.e. one-piece removal. Biatain CMC Alginate's fast gelling effect creates and maintains an optimal moist wound environment upon contact with exudate.

Latex free



### Biatain® Alginate

A dressing consisting of 85% calcium alginate and 15% carboxymethylcellulose (CMC). This structure is unique because the calcium alginate and CMC molecules are intertwined. This structure ensures that the wound exudate is absorbed vertically resulting in local absorption and hence minimising the risk of maceration.



| Code | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|---------|-------------------------|------------------|---------------------|
| 3705 | 29103705    | 5 x 5   | Sterile                 | 30               | 720                 |
| 3710 | 29103710    | 10 x 10 | Sterile                 | 10               | 240                 |
| 3715 | 29103715    | 15 x 15 | Sterile                 | 10               | 240                 |

### Biatain® Alginate (Filler)

Biatain Alginate filler is designed specifically for cavity wounds. The rope is presented pre-fluffed in order to maximise absorption capacity.



| Code | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|---------|-------------------------|------------------|---------------------|
| 3740 | 291037402   | 3 x 44  | Sterile                 | 6                | 120                 |

### Comfeel® Plus

Comfeel Plus dressings have a unique combination of Hydrocolloid and calcium alginate offering these benefits:



- Superior exudate management with fewer dressing changes and less disruption to wound
- Speciality shapes and sizes for tough to dress areas
- Colour guide indicates best time to change dressing
- 3-piece non-touch application system for ease of use
- · Latex free

### Comfeel® Plus

This innovative dressing incorporates the latest Hydrocolloid technology plus Calcium Alginate with a 'semi-permeable' top film to provide Superior Exudate Management via: better absorption, evaporation and adhesion. Features also include 'colour change indicator' and bevelled edges.



| Code  | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|-------|-------------|---------|-------------------------|------------------|---------------------|
| 33146 | 29103146    | 6 x 4   | Sterile                 | 30               | 720                 |
| 33110 | 29103110    | 10 x 10 | Sterile                 | 10               | 240                 |
| 33115 | 29103115    | 15 x 15 | Sterile                 | 5                | 240                 |
| 33285 | 29133285    | 17 x 17 | Sterile Sacral          | 5                | 80                  |
| 33120 | 29103120    | 20 x 20 | Sterile                 | 5                | 80                  |









### Comfeel® Plus Transparent

A thin, transparent Hydrocolloid with an intelligent semi-permeable top film. The high MVTR provides superior fluid handling capacity, resulting in longer wear time. A dressing that not only absorbs exudate but allows for visual inspection of the wound whilst in situ. The scar management quality makes this product the perfect choice for post-op wound management. Now comes with bevelled edges.



| Code  | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|-------|-------------|---------|-------------------------|------------------|---------------------|
| 33530 | 29103530    | 5 x 7   | Sterile                 | 10               | 300                 |
| 33548 | 29103548    | 5 x 25  | Sterile                 | 5                | 140                 |
| 33533 | 29103533    | 10 x 10 | Sterile                 | 10               | 240                 |
| 33536 | 29103536    | 9 x 14  | Sterile                 | 10               | 240                 |
| 33537 | 29103537    | 9 x 25  | Sterile                 | 5                | 140                 |
| 33539 | 29103539    | 15 x 15 | Sterile                 | 5                | 240                 |
| 33542 | 29103542    | 15 x 20 | Sterile                 | 5                | 120                 |
| 33545 | 29103545    | 20 x 20 | Sterile                 | 5                | 80                  |

### Comfeel® Plus Contour

This versatile dressing is designed for difficult to dress areas, combining a centre of Comfeel Plus dressing with a transparent Hydrocolloid framework.



| Code  | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|-------|-------------|---------|-------------------------|------------------|---------------------|
| 33280 | 29103280    | 6 x 8   | Sterile Heel            | 5                | 240                 |
| 33283 | 29103283    | 9 x 11  | Sterile Sacral          | 5                | 80                  |

### Comfeel® Plus Pressure Relief

Suitable for the prevention of pressure injuries in patients classified 'at risk'. The polyurethane foam pad on the back can be customised to suit the size of the area at risk.



| Code  | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|-------|-------------|---------|-------------------------|------------------|---------------------|
| 33350 | 6748        | 7       | Sterile Heel            | 10               | 180                 |
| 33353 | 6749        | 10      | Sterile Hip             | 10               | 180                 |

## Biatain® Ag

Major independent studies have conclusively proved that Biatain Ag is the only silver dressing that is effective on all bacteria commonly found in infected non healing wounds. Ensuring fast killing of all bacteria and effective exudate management, Biatain Ag is the best solution for infected wounds.



### Biatain® Ag provides:

- Broad spectrum silver profile proven to be effective against resistant bacteria
- Intelligent and sustained silver release
- Prevents wound infection
- Optimal wound healing environment
- Latex free

### Biatain® Ag Anti Bacterial Non-Adhesive Foam

Combines effective antimicrobial technology with superior absorption for faster healing of infected wounds.



| Code | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|---------|-------------------------|------------------|---------------------|
| 5105 | 13294       | 5 x 7   | Sterile                 | 5                | 160                 |
| 9622 | 11860       | 10 x 10 | Sterile                 | 5                | 240                 |
| 9625 | 13295       | 15 x 15 | Sterile                 | 5                | 80                  |

### Biatain® Ag Anti Bacterial Adhesive Foam

Additional benefit of a Hydrocolloid based adhesive border makes Biatain Aq adhesive an ideal dressing for infected wounds with healthy surrounding skin.



| Code | USL<br>Code | Size cm*    | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|-------------|-------------------------|------------------|---------------------|
| 9631 | 13296       | 7.5 x 7.5   | Sterile                 | 5                | 160                 |
| 9632 | 13297       | 12.5 x 12.5 | Sterile                 | 5                | 240                 |

<sup>\*</sup>Indicates dressing size, for pad size refer to retail box

### Biatain® Ag Cavity

Biatain Ag Cavity is indicated for exuding cavity wounds. Biatain Ag Cavity contains silver as an antibacterial agent for the treatment of critically colonised wounds and the revival of the healing process. Perforations in the dressing, allow for easy folding into a deep cavity. The dressing can also be easily cut to fit any wound.



| Code | USL<br>Code | Size cm | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|---------|-------------------------|------------------|---------------------|
| 9628 | 13301       | 5 x 8   | Sterile                 | 5                | 160                 |

### Purilon® Hydrogel

A balanced Hydrogel that provides fast and effective debridement. A unique match of hydrating and absorbing properties, combined with a high cohesive power, ensures that the Purilon Gel is the state of art Hydrogel.



| Code | USL<br>Code | Size     | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|----------|-------------------------|------------------|---------------------|
| 3906 | 29103906    | 8g Tube  | Sterile                 | 10               | 180                 |
| 3900 | 29103900    | 15g Tube | Sterile                 | 10               | 180                 |

### Comfeel® Wound Grid

| Code | USL<br>Code | Size cm   | Sterile/<br>Non Sterile | Units Per<br>Box | Units Per<br>Carton |
|------|-------------|-----------|-------------------------|------------------|---------------------|
| 3052 | 2 29103052  | 22 x 15.5 | Non Sterile             | 50               | 1000                |



FOR FURTHER INFORMATION CONTACT
USL Customer Services
FP 0800 658 814 FF 0800 830 660
e. customerservices@uslmedical.co.nz

Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare.

Our business includes ostomy care, continence care, wound and skin care and urology care. We operate globally and employ more than 10,000 employees.

Ostomy Care / Continence Care / Wound & Skin Care / Urology Care

